



Notes

1. Please read this manual carefully before doing experiment.
2. Please check whether there is crystallization or precipitation in the lysis buffer before use. If it happens, please dissolve it in a 56 °C water bath.

Viral DNA/RNA Kit

Product: Viral DNA/RNA Kit

Size: 50 tests/Box, 200 tests/Box

Applications

It is used for nucleic acids isolation, enrichment and purification. The processed products can be used for clinical diagnosis in vitro.

Principle

Nucleic acid isolation and purification kit provides a simple, rapid and efficient method to extract DNA/RNA from swab, serum, plasma, bronchoalveolar lavage fluid and other cell-free body fluids. The unique buffer system enables the nucleic acid in the lysate to be efficiently and specifically binded to the silica-based membrane. The obtained nucleic acid has high purity, stable quality, and is free of protein, nuclease and other contaminants and inhibitors. It can be applied to various conventional operations, including PCR, fluorescence quantitative PCR and other experiments.

Manufacturer

Jiangsu CoWin Biotech Co., Ltd
Add: No.58 G52 Building, East of Lujia Road, West Side of Tai Road, China medical city, Taizhou, Jiangsu, 225300
Phone: 86-523-86201352 Fax: 86-523-86816890
Web: <https://www.cwbiotech.com> E-mail: service@cwbiotech.com

EC-Representative

SUNGO Europe B.V.
VAT: NL857821659B01
Add: Olympisch Stadion 24, 1076DE Amsterdam, Netherlands
Contact: SUNGO Secretary
Tel /Fax: +31 (0) 2021 11106
E-mail: ec.rep@sungogroup.com

Version: A/2

Effective Date: 2020-06-05

Kit Components

Component	50 tests/Box		200 tests/Box	
	Size	Quantity	Size	Quantity
Lysis Buffer	15 ml/bottle	1	50 ml/bottle	1
Washing Buffer 1	30 ml/bottle	1	120 ml/bottle	1
Washing Buffer 2	30 ml/bottle	1	120 ml/bottle	1
RNase-Free Water	10 ml/bottle	1	25 ml/bottle	1
Adsorption Columns	50/package	1	50/package	4
Collection Tubes	50/package	1	50/package	4

Storage Condition and Valid Period

Store at 0-35°C up to 12 months.

We suggest products transportation at 0-40°C for no more than 7 days.

Sample Requirements

Applicable Samples: swab, serum, plasma, bronchoalveolar lavage fluid and other cell-free body fluids.

Procedure

Equipment and Reagents to Be Supplied by User: Constant temperature mixer (CW BIO, CW2593 is recommended), isopropanol.

Things to do before starting: Mix all the reagents and gently invert 3-5 times before use.

General purification procedure

1. Take 1.5 ml centrifugal tube (supplied by user), add 200 µl sample (the sample needs to be balanced to room temperature), 200 µl lysis buffer, 300 µl isopropanol. Vortex for 5 seconds, and store it at room temperature and vortex for 10 minutes with a constant temperature mixer of 1200 rpm.

Note: For wet swab, take 200 µl sample, mix it thoroughly. Soak the dry swab sample in 400 µl normal saline and mix well, let stand for 5 minutes, centrifugate at 12000 rpm for 1 minute, and extract 200 µl.

2. The obtained solution from step 1 was transferred to the adsorption column of collection tube. Centrifugation at 12000 rpm (~ 13400 g) for 1 minute, and discard the liquid. And then put the adsorption column back into the collection tube.
3. Add 500 µl washing buffer 1 to the adsorption column, centrifugate at 12000 rpm for 1 minute, discard the liquid, and put the adsorption column back into the collection tube.
4. Add 500 µl washing buffer 2 to the adsorption column, centrifugate at 12000 rpm for 1 minute, discard the liquid, and put the adsorption column back into the collection tube.
5. Centrifugation at 12000 rpm for 2 minutes, and discard waste liquid in the collection tube. Place the column at room temperature for 2 minutes and let it dry.
6. Add 40-100 µl RNase-Free water to the middle part of the adsorption column membrane, keep 2 minutes at room temperature and centrifugate for 1 minute at 12000 rpm, the RNA solution was preserved at -80°C to prevent degradation.

Rapid purification procedure (swab sample only)

1. Take 1.5 ml centrifugal tube (supplied by user), add 300 µl sample (the sample needs to be balanced to room temperature), 200 µl lysis buffer, 300 µl isopropanol. Vortex for 5 seconds, and centrifugate for a second to ensure no liquid residue on the lid and wall of the tube.
Note: For wet swab, take 300 µl sample, mix it thoroughly. For dry swab, soak the sample in 500 µl normal saline and mix well, let stand for 5 minutes, centrifugate at 12000 rpm for 1 minute, and extract 300 µl.
2. The obtained solution from step 1 was transferred to the adsorption column. Centrifugation at 12000 rpm (~ 13400 g) for 30 seconds, and discard the liquid. And then put the adsorption column back into the collection tube.
3. Add 500 µl washing buffer 2 to the adsorption column, centrifugate at 12000 rpm for 30 seconds, discard the liquid, and put the adsorption column into a new collection tube.
4. Add 40-100 µl RNase-Free water to the middle part of the adsorption column membrane, centrifugate for 30 seconds at 12000 rpm, the RNA solution was preserved at -80°C to prevent degradation.